



ExHC

Excel Formula®

Heat Cure Denture Base Material

Finally, a higher impact heat-cure denture base that combines superior workability with a quick cure option! Excel Formula® Heat Cure denture base is a lab technician's dream. Guaranteed quality performance at competitive rates. Saves you time, effort, and money.

- ◆ 3 curing options: quick cure, overnight or microwavable
- ◆ Great impact strength
- ◆ Reaches snap dough consistency within 5 minutes, workable for another 30min.
- ◆ Highly color stable with matching repair and relining resins
- ◆ Super-fine particle size enhances strength and eliminates porosity
- ◆ Low residual monomer



Light Veined Veined Reddish Veined Dark Veined Clear

Available Products:



Universal Pack - 1 Kilo (2.2 lbs) powder, 500 mL (1 pint) liquid and directions for use

- 10111-Veined ■ 10112-Light Veined ■ 10113-Reddish Veined ■ 10114-Clear ■ 10115-Dark Veined



1 Kilo Powder Only (2.2 lbs.)

- 10211-Veined ■ 10212-Light Veined ■ 10213-Reddish Veined ■ 10214-Clear ■ 10215-Dark Veined



10 Kilo Powder Only (22 lbs.)

- 10221-Veined ■ 10222-Light Veined ■ 10223-Reddish Veined ■ 10224-Clear ■ 10225-Dark Veined



5 x 10 Kilo Powder Only (110 lbs.)

- 10231-Veined ■ 10232-Light Veined ■ 10233-Reddish Veined ■ 10234-Clear ■ 10235-Dark Veined



500 mL Liter Liquid Only

- 10310



1 Liter Liquid Only

- 10340



4 x 1 Liter Liquid Only

- 10350



5 Liter Liquid Only

- 10320

