

ExHIP

Excel Formula[®] Hi-Impact Pourable

World's first dental acrylic to achieve high strength and high impact resistance from a cold cure formulation. More than 4x stronger and 60% tougher than conventional cold-cure denture base polymers. Offers a 75% time reduction from mixing to removal from flask, improving lab productivity and reduces costs.

- ◆ High impact strength
- ◆ Exceptional flexural strength
- ◆ 4-5 minute pour time/20 minute short curing time
- ◆ Super-fine particle size enhances strength and tooth bond
- ◆ Cadmium free
- ◆ Easy to mix to a perfect consistency
- ◆ Molecular structure eliminates porosity and minimizes shrinkage
- ◆ Exceptional color stability
- ◆ Silicon free



Available Products:



Universal Pack - 1 Kilo (2.2 lbs.) Powder, 500 mL (1 Pint) Liquid

■ 24111-Veined 99



1 Kilo Powder Only (2.2 lbs.)

■ 24211-Veined 99



10 Kilo Powder Only (22 lbs.)

■ 24221-Veined 99



500ml Liter Liquid Only

■ 24330




1 Liter Liquid Only

■ 24340



Veined 99

FRACTURE TOUGHNESS



FRACTURE WORK

