



ExP

Excel Formula[®] Lab Putty Dental Material

Excel Formula[®] Lab Putty introduces yet another level of performance and excellence to the lab. With a great chemical makeup and Shore hardness of 85, you are guaranteed to receive strong and accurate results time after time.

Available Products:



1.6 Kilo Base Putty (3.52 lbs.) and Scoop

■ 30101



10 Kilo Base Putty (22 lbs.)

■ 30120



60mL Catalyst

■ 30140



Pink

- ◆ Fine-grained material providing great detail and accuracy
- ◆ Offers a shore hardness of 85
- ◆ Catalyst blends easily into putty leaving hands grease-free
- ◆ 2-3 minutes of working time / sets in 8-10 minutes
- ◆ Premium material at a competitive price

