



Exp

Excel Formula® Pourable Denture Base Material

State-of-the-science polymer now brings unprecedented performance to pourable acrylics. With the Excel Formula® Pourable material, lab technicians can at last enjoy a pourable denture base that mixes quickly to produce the perfect consistency needed for optimum flow and unparalleled performance.

- ◆ 4-5 minute pour time/20 minute short curing time
- ◆ Super-fine particle size enhances strength and tooth bond
- ◆ Cadmium free
- ◆ Easy to mix to a perfect consistency
- ◆ Molecular structure eliminates porosity and minimizes shrinkage
- ◆ Exceptional color stability



Light Veined Veined Reddish Veined Dark Veined Clear

Available Products:



Universal Pack - 1 Kilo (2.2 lbs.) Powder, 500 mL (1 Pint) Liquid

- 23111-Veined ■ 23112-Light Veined ■ 23113-Reddish Veined ■ 23114 Clear ■ 23115 Dark Veined



1 Kilo Powder Only (2.2 lbs.)

- 23211-Veined ■ 23212-Light Veined ■ 23213-Reddish Veined ■ 23214 Clear ■ 23215 Dark Veined



10 Kilo Powder Only (22 lbs.)

- 23221-Veined ■ 23222-Light Veined ■ 23223-Reddish Veined ■ 23224 Clear ■ 23225 Dark Veined



1 Liter Liquid Only

- 23340



5x1 Liter Liquid Only

- 23350

